## SYRINGE SHIELD FOR SAFETY SYRINGES

Handle nuclides and help prevent needle sticks without compromising efficiency or radiation safety.



- Allows full function of safety syringes needle can be enclosed in the protective sleeve while syringe is in the shield
- Tungsten shielding with lead glass window
- Push-button release for one-handed syringe disposal
- Accommodates the following syringes:
   B-D Safety-Lok<sup>™</sup>
   Monoject<sup>®</sup> Safety

The Pro-Tec® V is designed to comply with OSHA needle-stick prevention guidelines without compromising radiation protection or procedural efficiency. This unique design allows the syringe's protective sleeve to be employed while the syringe is in the shield. It also provides a simple one-handed method to safely release the syringe from the shield. The barrel of the shield is constructed of .08" (2 mm) thick tungsten that will reduce radiation exposure from

thick tungsten that will reduce radiation exposure from Tc-99m by 94%. A lead glass window provides a shielded view of the injection fluid and the syringe markings. A white reflective surface on the shield interior improves viewing of the syringe's markings and fluid content. The window is recessed in its beveled mount to reduce scratching and prevent breakage.



## Here's how it works:

- Prior to injection: The safety syringe is pushed into the shield until it engages. At this point, a springloaded catch engages the syringe and holds it in place.
- 2. After injection: While the safety syringe is still in the shield, the syringe's protective sleeve is positioned over the needle and locked in place. Finger pressure on the "LOCK" side of the shield's dualaction button prevents the syringe's outer sleeve from turning within the shield. This allows the syringe flanges to be rotated to lock the protective sleeve.
- 3. Syringe disposal: While holding the shield over a sharps container, the "RELEASE" side of the button is pushed to drop the syringe directly into the container.

## SPECIFICATIONS Shielding:

**Barrel:** 0.08" thick (2mm) tungsten **Lead Glass Window:** 0.17" thick (4.3mm) 5.05 density (1.5mm lead equivalency)

## Pro-Tec<sup>®</sup> V Syringe Shield for Safety Syringes: Accommodates B-D Safety-Lok and Monoject Safety Syringes

<b>007-737</b> Syringe Shield, 3 cc	\$200.00
007-743 Syringe Shield, 5 cc	210.00
007-744 Syringe Shield, 10 cc	260.00

Prices subject to change without notice.

www.jgravengaard.com 2005-12-01