

# SHARPS CONTAINER SHIELDS

*Key-lock design meets OSHA standards*



Sharps Container Shields are a simple way to dispose of used syringes that may contain low energy gamma radiation residue. The shields are constructed of steel lined with .125" (.32 cm) lead. They feature a top with sliding port and a swing handle that allows easy transport.

Two models are available for one sharps container:

- Model 039-325 accommodates small and medium Monoject Sharps Containers.
- Model 039-335 accommodates Winfield and B-D Guardian Sharps Containers.

One model is available for two sharps containers:

- Model 039-350 accommodates two small Monoject Sharps Containers simultaneously for decay rotation.

The sharps containers are stacked inside the shield. When the top container is filled, the bottom decayed container is removed for proper waste disposal. The top container is then shifted to the lower position for decay while a fresh container is placed on top for immediate use. Spring clamps hold containers securely in place.



*Prices subject to change without notice.*

## SPECIFICATIONS:

*039-325 Sharps Shield. Monoject Shield uses small and medium Monoject Sharps containers*

**Dimensions:** 7.5" l x 11.75" w x 10" h (19 x 30 x 25 cm)

**I.D.:** 7" l x 11.25" w x 9.5" h (18 x 29 x 24 cm)

**Lead Shielding:** .125" thick (3 mm)

**Security:** Key locked

**Weight:** 48 lb (22 kg)

*039-335 Sharps Shield. Winfield Shield uses Winfield and B-D Guardian Sharps containers*

**Dimensions:** 7.5" l x 12.5" w x 13" h (19 x 32 x 33 cm)

**I.D.:** 6.875" l x 10.75" w x 12.5" h (17.5 x 27 x 32 cm)

**Lead Shielding:** .125" thick (3 mm)

**Security:** Key locked

**Weight:** 53 lb (24 kg)

*039-350 Sharps Shield. Dual Monoject Shield uses two small Monoject Sharps containers*

**Dimensions:** 7.375" l x 12" w x 16.75" h (19 x 30.5 x 42.5cm)

**I.D.:** 7" l x 11.5" w x 14.875" h (18 x 29 x 38 cm)

**Lead Shielding:** .125" thick (3 mm)

**Weight:** 60 lb (27 kg)



## Sharps Container Shields:

**039-325** Sharps Container Shield, .125" lead..... \$510.00  
*For Monoject containers 039-338 and 039-387.*

**039-335** Sharps Container Shield, .125" lead..... 510.00  
*For Winfield container 039-341 and B-D Guardian container 039-321.*

**039-350** Sharps Shield Container, Dual, .125" lead... 650.00  
*For Monoject containers 039-338.*

## Related:

**039-338** Monoject Container, small, 10/pkg ..... \$50.00  
*Fits 039-325. 039-350 and 039-330.*

**039-387** Monoject Container, medium, 10/pkg ..... 85.00  
*Fits 039-325.*

**039-341** Winfield Container, medium (#182), 12/pkg... 75.00  
*Fits 039-335.*

**039-321** B-D Guardian Container, 8/pkg ..... 85.00  
*Fits 039-335.*

[www.jgravengaard.com](http://www.jgravengaard.com) 2007-07-09

**J. Gravengaard Corporation**

Medical Imaging Sales & Service

PO Box 20171 Portland, Oregon 97294 Nationwide (800)796-SCAN Local (503)255-0535 Fax (503)255-7047